

Swept Source Optical Coherence Tomography



Swept Source OCT Imaging System shown on a Thorlabs ScienceDesk™. See page 26 for details on Thorlabs' ScienceDesks.

SS-OCT Features

- 1040 or 1325 nm Center Wavelength
- Ideal for Imaging Small Animals, Tissue, and Other Highly Scattering Material
- Video Microscopy Imaging with Integrated CCD Camera in Handheld Probe
- Doppler Flow Imaging Included
- Compact System Design Ready for OEM Applications
- 3D Imaging Capability
- Sub-Surface Imaging up to 3 mm in Tissue

Optional Features

- Family of Objective Lenses to Optimize Imaging Performance
- Polarization-Sensitive OCT
- High-Power Light Sources
- Special Center Wavelength Sources for Industrial Samples



Swept Source Optical Coherence Tomography (SS-OCT) is a Fourier Domain OCT technique that utilizes a frequency-swept, broadband tunable laser and balanced detection scheme to generate and measure frequency-based interferometric signals. With SS-OCT the spectral components of the interferogram is encoded in time as the laser source tunes across a broad bandwidth. Compared to Spectral Domain OCT technology, SS-OCT systems do not suffer from inherent sensitivity degradation at greater depths. Therefore, SS-OCT is preferred for applications where deeper imaging is desired.

Biological Imaging

Three-Dimensional Imaging of Frog Tadpole

The data shown here represents a series of *en face* images of an African frog tadpole and a 3D image reconstruction. The data was acquired using the Thorlabs OCS1300SS OCT microscope system. All images are 6 mm x 8 mm and were taken from the posterior to the anterior of the tadpole in 100 μ m increments.

Reference: R. Huber, *et al.*, *Optics Express* **13**, 10523 (2005)



Swept Source OCT Systems

OCS1300SS: Deep Imaging

The OCS1300SS is our most popular OCT imaging system. Utilizing a 55 kHz swept laser source with a center wavelength of 1325 nm (SL1325-P55), this system provides deep imaging (up to 3 mm) in highly scattering tissue. The OCS1300SS comes with three-dimensional, as well as Doppler imaging capabilities bundled into a complete data control and acquisition software package. The OCS1300SS Imaging System includes a 3D imaging capable handheld probe with articulating stand. Imaging performance is enhanced by using an OCT-optimized objective lens (Thorlabs' Item# LSM03). Details on this and other OCT objectives can be found on page 1386.

OCS1050SS: High-Resolution Imaging

The OCS1050SS builds upon the standards set by our 1300 nm system by offering high-resolution swept-source OCT images while utilizing a wavelength band that is less susceptible to water absorption. It is capable of imaging depths up to 2.5 mm, with an axial resolution of 12.5 μm (in air). The system includes the same 3D imaging handheld probe; when coupled with an LSM02-BB scan lens, the OCS1050SS is Thorlabs' highest resolution swept source OCT system. Our new 1050 nm swept laser source (SL1050-P55) has an 80 nm bandwidth, 5 mm coherence length, and 10 mW output power. See page 1368 for details on our high-speed, 55 kHz swept laser sources.

SS-OCT Systems	1300 nm	1050 nm
	■ Deep Imaging ■ A-Scan Rate: 55 kHz ■ Image Acquisition Rate: 50 Frames per Second**	■ High Resolution ■ A-Scan Rate: 55 kHz ■ Image Acquisition Rate: 50 Frames per Second**
Center Wavelength	1325 $\pm$ 15 nm	1040 $\pm$ 15 nm
Spectral Bandwidth (20 dB)	120 $\pm$ 10 nm	80 $\pm$ 10 nm
Axial Resolution (Air/Water)	12 μm / 9.0 μm	10 μm / 7.5 μm
Lateral Resolution	15 μm	12.5 μm
Maximum Imaging Depth*	6.0 mm	5.0 mm
Maximum Pixels per A-Scan	1024	1024
Optical Power on Sample (Average)	5.0 mW	4.0 mW


*Maximum Imaging Depth that is supported by the data acquisition system and software. 6.0 mm (5.0 mm for 1050 SS-OCT) imaging depth is susceptible to signal degradation due to finite coherence length of the swept laser source. Actual Imaging Depth will vary based on absorption and scattering properties of the sample. **Image Acquisition Rate based on 1024 lines per frame.

Brain Functional Imaging

The pictures to the right show OCT images representing the functionally active region of a rat brain. Images were taken by researchers at MIT and Massachusetts General Hospital using a modified Thorlabs Swept Source OCT system.

The OCT scan is directed to the region of interest using a video microscope. The OCT image enables identification of the skull (S), surface vasculature (V), and meningeal layers, including the dura mater (D). The blue-colored areas show the functionally active regions of the cortex, which were co-registered on the CCD camera and OCT channels.

Reference: A.D. Aguirre Y. Chen, J.G. Fujimoto, L. Ruvinskaya, A. Devor, and D.A. Boas, "Depth-resolved imaging of functional activation in the rat cerebral cortex using optical coherence tomography," *Optics Letters* **31**, 3459 (2006).



Doppler OCT Imaging

Doppler OCT is an extension of OCT that enables imaging of particle motion within a sample. In Fourier Domain OCT (FD-OCT) systems, there are no additional hardware requirements for implementation of Doppler imaging. Doppler OCT imaging capability is embedded in the software provided with all Thorlabs' SS-OCT systems and is ideal for studying embryonic cardiac dynamics, functional vascular imaging, or monitoring vascular treatment response. It is also useful for general flow velocimetry used in microfluidic channel monitoring.

Principles of Doppler OCT

The Fourier transform of the interferogram acquired in FD-OCT imaging (A-scan) produces a complex signal $[I(z) + iQ(z)]$, where the magnitude of that signal is used to create the structural OCT image. The complex portion of the signal contains information based on the phase of the interferogram. Any change in phase between consecutively acquired A-scans can be attributed to a Doppler frequency shift induced by particle motion. In Thorlabs' implementation of Doppler OCT, Doppler frequency shifts are calculated based on spatially averaging phase shifts within a sliding 2D window and

using a Kasai autocorrelation function (see Ref. 1). The Doppler frequency shift f_D caused by moving particles is related to the phase shift between A-scans, as described in the following expression:

$$f_D = \frac{\Delta\phi}{2\pi} f_A$$


where $\Delta\phi$ is the average phase shift within the sliding 2D window and f_A is the A-scan rate of the OCT system. The Thorlabs Doppler Imaging Mode displays the phase shift induced by moving particles using a standard Doppler colormap where red to yellow (violet to blue) indicates flow in the direction toward (away from) the OCT sample beam.

The mean velocity $\langle v \rangle$ of the moving particles, at any depth, can be quantified by knowing θ , the angle between the OCT sample beam and flow vector:

$$\langle v \rangle = \frac{\lambda_0 f_D}{2n \cos(\theta)}$$



Here, λ_0 is the center optical wavelength of the OCT sample beam and n is the index of refraction of the sample.

[1] C. Kasai, K. Namekawa, A. Koyano, R. Omoto *et al.* "Real-time two-dimensional blood flow imaging using an autocorrelation technique," *IEEE Trans. Sonics. Ultrason.* 32 458-464 (1985).



Developmental Biology

Thorlabs' Swept Source OCT Imaging System (OCS1300SS) with Doppler Imaging was used by researchers at the University of Toronto to study the cardiovascular system of living tadpoles. The series of images below show *in vivo* cross-sectional SS-OCT images of a beating tadpole heart superimposed with Doppler blood flow images. An optical Doppler cardiogram was obtained using a gated technique to increase the effective frame rate and improve the signal-to-noise ratio. The gating technique provides ultra high-speed visualization of the heart blood flow pattern in developing African frog embryos in both 3D and 4D (i.e., 3D + time) modes. This allows detailed visualization of the complex cardiac motion and hemodynamics in the beating heart.



Reference: A. Mariampillai, B.A. Standish, N.R. Munce, C. Randall, G. Liu, J.Y. Jiang, A.E. Cable, I.A. Vitkin, V.X.D. Yang, "Doppler optical cardiogram gated 2D color flow imaging at 1000 fps and 4D *in vivo* visualization of embryonic heart at 45 fps on a swept source OCT system", *Optics Express* 15, 1627 (2007).





Fig. A shows the 3D surface reconstruction of the tadpole heart, while Fig. B demonstrates the complex blood flow pattern of the heart via a 3D color Doppler map.

Swept Source OCT System Specifications

NEW
product

Due to rapid improvements to our systems, all specifications presented here are subject to change. Please refer to our website, www.thorlabs.com, for or call for the latest specifications



Handheld Probe (above) and
Articulated Probe Mount
(below) Included with all
OCT Systems



Optical*	OCS1300SS	OCS1050SS
Center Wavelength	1325 nm	1040 nm
Spectral Bandwidth (20 dB)	120 nm	80 nm
A-Scan Rate	55 kHz	55 kHz
Average Optical Power at Sample	5.0 mW	4.0 mW
Electronics		
A/D Conversion Rate	180 MS/s	
A/D Resolution	16 Bit	
A/D Channels	2	
Analog Output Rate	1 MS/s	
Analog Output Resolution	16 Bit	
Analog Output Channels	4	
Imaging Specifications		
2D Cross-Sectional OCT Imaging		
Imaging Speed (1024 A-Scans/Frame)	50 frames per second	
Maximum Sampling Resolution	8192 (W) x 1024 (D) Pixels	
Maximum Scan Range	10 mm	
Maximum Imaging Depth**	6.0 mm	5.0 mm
Lateral Resolution with LSM03 Objective	15 μm	12.5 μm
Axial Resolution in Air (Water)	12 μm (9.0 μm)	10 μm (7.5 μm)
3D Volume OCT Imaging		
Maximum Volume Size**	10 mm x 10 mm x 6.0 mm	10 mm x 10 mm x 5.0 mm
Maximum Sampling Resolution	1024 (W) x 1024 (L) x 1024 (D)	
Volume Acquisition Time	5 s (512 x 512 x 1024)	
En Face Video Microscopy		
CCD Camera	0.25 Megapixels (Effective)	
Maximum Pixel Resolution	328 x 256	
Imaging Speed	29 Frames/s	

*Longer coherence length system configurations are available with shorter turning ranges, please call for details.

**Imaging Depth is software selectable up to the maximum range.

ITEM#	\$	£	€	RMB	DESCRIPTION
OCS1300SS	\$ 60,000.00	£ 41,594.00	€ 53,269.00	¥ 506,642.00	1300 nm SS-OCT Imaging System
OCS1050SS	\$ 60,000.00	£ 41,594.00	€ 53,269.00	¥ 506,642.00	1050 nm SS-OCT Imaging System