

56 Sparta Avenue • Newton, New Jersey 07860
 (973) 300-3000 Sales • (973) 300-3600 Fax
 www.thorlabs.com

THORLABS

M54S76 - December 19, 2014

Item # M54S76 was discontinued on December 19, 2014. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

FIBER OPTIC TERMINATION TOOLS

- ▶ Buffer Stripping Tools for Ø80 to Ø1550 µm Clad Fibers
- ▶ Furcation Stripping Tools for Ø900 µm to Ø3.8 mm Fiber Jackets

T06S13
Fiber Stripping Tool



FTS3
Stripping Tool for Ø3 mm and
Ø3.8 mm Furcation Tubing



CT042
Connector
Crimp Tool



FTS4
3 Hole Fiber Stripper



AFS900
Adjustable Fiber and
Wire Stripping Tool



T865
Kevlar Cutting
Shears



[Hide Fiber Buffer Stripping Tools](#)


Fiber Buffer Stripping Tools

[Hide Overview](#)

OVERVIEW

Features

- Self-Aligning Blade Set Assures Concentric Scoring of Buffer or Coating
- Color-Coded Blades are Long-Lasting and Swappable
- Each Stripping Tool Includes a BFG1 Bare Fiber Gripper
- Foolproof, No-Nick Design
- Fast, Reliable Fiber Stripping



Click to Enlarge
 T12S16 Stripping Tool shown with BFG1 Fiber Gripper

These fiber buffer stripping tools provide a quick, easy, and reliable way to remove the buffer from an optical fiber in preparation for connectorization. A fiber guide and matched blades ensure that the optical fiber is correctly positioned and stripped each time. The blades are color coded to allow for fast identification of the proper fiber stripping tool. All fiber stripping tools include the BFG1 bare fiber gripper.

Fiber Gripper

The BFG1 bare fiber gripper is an ideal aid when removing buffer from an optical fiber. The textured rubber material of this gripper provides a secure hold without damaging the fiber. This method of gripping is recommended over wrapping or clamping the loose end of the fiber, techniques which can create

microfractures in the fiber.

[Hide Specs](#)

S P E C S

Tool Selection

Step 1: Note your cladding and coating diameters along with their respective \pm tolerances.

Step 2: Look down the second column of the table below for your fiber size.

Step 3: With your fiber size identified in the chart below, scan across to the corresponding 'Cladding Range' and 'Coating Range' columns. Ensure that your fiber dimensions \pm tolerances fall within the range listed. If they do, then refer to the corresponding Item # to place your order. If the maximum fiber dimensions fall outside the range shown, go to the next larger tool.

Item #	Typical Fiber Cladding/Coating	Cladding Range	Coating Range ^a
T04S10	80 μ m / 170 μ m	65 - 80 μ m	150 - 250 μ m
T06S13	125 μ m / 250 μ m	125 - 135 μ m	250 - 343 μ m
T06S16	125 μ m / 400 μ m	125 - 135 μ m	343 - 407 μ m
T08S13	125 μ m / 250 μ m	125 - 175 μ m	250 - 343 μ m
T08S40	125 μ m / 900 μ m	125 - 175 μ m	889 - 1016 μ m
T10S13	200 μ m / 300 μ m	180 - 230 μ m	250 - 343 μ m
T12S16	230 μ m / 400 μ m	235 - 280 μ m	343 - 407 μ m
T12S18	230 μ m / 430 μ m	235 - 280 μ m	407 - 457 μ m
T12S21	230 μ m / 500 μ m	235 - 280 μ m	457 - 533 μ m
T12S25	230 μ m / 600 μ m	235 - 280 μ m	533 - 635 μ m
T14S21 ^b	285 μ m / 500 μ m	285 - 330 μ m	457 - 533 μ m
T16S31	325 μ m / 650 μ m	335 - 380 μ m	635 - 787 μ m
T18S25	400 μ m / 600 μ m	385 - 430 μ m	533 - 635 μ m
T18S31	400 μ m / 730 μ m	385 - 430 μ m	635 - 787 μ m
T21S31	425 μ m / 730 μ m	435 - 500 μ m	635 - 787 μ m

T23S46	500 μm / 1000 μm	505 - 550 μm	1016 - 1168 μm
T28S46	630 μm / 1040 μm	605 - 680 μm	1016 - 1168 μm
M37S46	830 μm / 1040 μm	835 - 900 μm	1016 - 1168 μm
M37S63	830 μm / 1400 μm	835 - 900 μm	1397 - 1600 μm
M44S63	1035 μm / 1400 μm	905 - 1050 μm	1397 - 1600 μm
M44S67	1035 μm / 1600 μm	905 - 1050 μm	1600 - 1702 μm
M54S76	1240 μm / 1650 μm	1055 - 1350 μm	1778 - 1930 μm
M63S86	1550 μm / 2000 μm	1390 - 1600 μm	2057 - 2184 μm

- Coating refers to the jacket, buffer, or coating that is being removed.
- T14S21 is the only strip tool compatible with BFL48-200 and BFH48-200 fiber.

Part Number	Description	Price	Availability
T04S10	Fiber Stripping Tool, Typical Cladding/Coating: 80 μm / 170 μm	\$135.66	Today
T06S13	Fiber Stripping Tool, Typical Cladding/Coating: 125 μm / 250 μm	\$69.36	Today
T06S16	Fiber Stripping Tool, Typical Cladding/Coating: 125 μm / 400 μm	\$69.36	Today
T08S13	Fiber Stripping Tool, Typical Cladding/Coating: 125 μm / 250 μm	\$69.36	Today
T08S40	Fiber Stripping Tool, Typical Cladding/Coating: 125 μm / 900 μm	\$69.36	Today
T10S13	Fiber Stripping Tool, Typical Cladding/Coating: 200 μm / 300 μm	\$69.36	Today
T12S16	Fiber Stripping Tool, Typical Cladding/Coating: 230 μm / 400 μm	\$66.61	Today
T12S18	Fiber Stripping Tool, Typical Cladding/Coating: 230 μm / 430 μm	\$66.61	Today
T12S21	Fiber Stripping Tool, Typical Cladding/Coating: 230 μm / 500 μm	\$66.61	Today
T12S25	Fiber Stripping Tool, Typical Cladding/Coating: 230 μm / 600 μm	\$66.61	Today
T14S21	Fiber Stripping Tool, Typical Cladding/Coating: 285 μm / 500 μm	\$66.61	Today
T16S31	Fiber Stripping Tool, Typical Cladding/Coating: 325 μm / 650 μm	\$65.28	Today
T18S25	Fiber Stripping Tool, Typical Cladding/Coating: 400 μm / 600 μm	\$66.61	Today
T18S31	Fiber Stripping Tool, Typical Cladding/Coating: 400 μm / 730 μm	\$65.28	Today
T21S31	Fiber Stripping Tool, Typical Cladding/Coating: 425 μm / 730 μm	\$65.28	Today
T23S46	Fiber Stripping Tool, Typical Cladding/Coating: 500 μm / 1000 μm	\$65.28	Today
T28S46	Fiber Stripping Tool, Typical Cladding/Coating: 630 μm / 1040 μm	\$67.22	Today
M37S46	Fiber Stripping Tool, Typical Cladding/Coating: 830 μm / 1040 μm	\$67.83	Today
M37S63	Fiber Stripping Tool, Typical Cladding/Coating: 830 μm / 1400 μm	\$67.83	Today
M44S63	Fiber Stripping Tool, Typical Cladding/Coating: 1035 μm / 1400 μm	\$67.83	Today
M44S67	Fiber Stripping Tool, Typical Cladding/Coating: 1035 μm / 1600 μm	\$67.83	Today
M54S76	Fiber Stripping Tool, Typical Cladding/Coating: 1240 μm / 1650 μm	\$67.83	Lead Time
M63S86	Fiber Stripping Tool, Typical Cladding/Coating: 1550 μm / 2000 μm	\$79.25	Today
BFG1	Rubber Gripper for Bare Fiber Optics	\$5.00	Today

[Hide Crimping Tool for Connectors](#)

Crimping Tool for Connectors

[Hide Overview](#)

OVERVIEW

Fiber Connector Crimp Tool

This fiber connector crimp tool can be used for crimping SMA, FC, SC, and ST connectors. The crimp tool featured here has hex crimp sizes of: 0.128", 0.151", 0.178", 0.197", and 0.215".

Connectors with Ø3 mm or greater tubing require the crimp tool. The Ø900 µm tubing or smaller does not need to be crimped.

Part Number	Description	Price	Availability
CT042	Crimp Tool Handle with Die for SMA, FC, or ST Connector	\$99.00	Today

[Hide Stripping Tools](#)

Stripping Tools

[Hide Overview](#)

OVERVIEW

Kevlar Cutting Shears

The T865 cutting shears have carbon steel serrated blades that are specifically designed to cut through Kevlar fibers. The serrated blade keeps the Kevlar fibers from sliding along the blades as the shears are being closed, which allows the Kevlar fibers to be cut easily. It is important to only use the T865 cutting shears for cutting Kevlar, any other material, including the silica fiber itself, can potentially damage to the precision edge of the blades and make them less effective in cutting Kevlar. These shears feature padded handles that have a return spring.

Adjustable Fiber Buffer and Jacket Stripping Tool

The AFS900 stripping tool is a versatile tool capable of stripping both the fiber jacket (furcation tubing) and fiber buffer. The blades have a V-groove slot to precisely hold the fiber in the proper position as the stripping tool is closed. The AFS900 has an adjustable blade stop that can be used to ensure that the optical fiber is not cut when stripping the jacket or buffer. This tool is often used when stripping $\text{\O}900\ \mu\text{m}$ fiber jacket (tight or loose) or a $\text{\O}250\ \mu\text{m}$ buffer from a $\text{\O}125\ \mu\text{m}$ clad fiber. Because of the delicate nature of fiber, we recommend that when the AFS900 is used to strip fiber buffer that the stop be properly set and then to dedicate the tool for that purpose. Also available (see above) are fiber buffer stripping tools with blades that are dedicated for use with specific fiber cladding and buffer diameters.

Stripping Tool for $\text{\O}3\ \text{mm}$ and $\text{\O}3.8\ \text{mm}$ Furcation Tubing

The FTS3 will provide precise cuts through $\text{\O}3\ \text{mm}$ and $\text{\O}3.8\ \text{mm}$ furcation tubing. Simply set the depth of the dual cutting blades before placing the cutting tool around the fiber. Then spin the tool around the fiber several times. This results in two precise cuts through the tubing. The furcation tubing between the two cuts can then be removed by using the blade at the end of the FTS3 to make a precise cut through the tubing along the length of the fiber. Each FTS3 also comes with a blade that has a rounded tip. Long cuts along the length of the fiber can be made through the furcation tubing by swapping the blade at the end of the FTS3 with the rounded blade. Simply place the fiber in the end V-groove, press the blade through the furcation tubing, and then draw the fiber through the tool. Replacement blades (FTS3B) are available for the FTS3 stripping tool. Each FTS3B contains three straight blades plus the round tipped blade.

Three-Hole Fiber Stripping Tool

The FTS4 fiber stripping tool performs all common fiber stripping functions. The first hole, indicated by the 1 in the photo to the right, strips the $\text{\O}1.6 - 3\ \text{mm}$ fiber jacket down to the $\text{\O}600 - 900\ \mu\text{m}$ buffer coating. The second hole strips the $\text{\O}600 - 900\ \mu\text{m}$ buffer coating down to the $\text{\O}250\ \mu\text{m}$ coating. The third hole is used to strip the $\text{\O}250\ \mu\text{m}$ cable down to the glass fiber without nicks or scratches.

The handle is made of TPR (Thermoplastic Rubber) and the tool is 6" long from the tip of the blade to the end of the handle. This fiber stripper is preset at the factory and cannot be adjusted or calibrated. The stripping tool should be periodically checked for proper operation, and can be cleaned with a fine plastic bristled brush to remove debris from the openings. The openings can also be cleaned with dry compressed air CA3.



Click to Enlarge

Fiber Gripper

The BFG1 bare fiber gripper is an ideal aid when removing buffer from an optical fiber. The textured rubber material of this gripper provides a secure hold without damaging the fiber. This method of gripping is recommended over wrapping or clamping the loose end of the fiber, techniques which can create microfractures in the fiber.

Part Number	Description	Price	Availability
AFS900	Adjustable Stripping Tool for Fiber Buffers, Fiber Jackets, and Wire	\$16.90	Today
T865	Kevlar Cutting Shears	\$36.67	Today
FTS3	Stripping Tool for $\text{\O}3\ \text{mm}$ and $\text{\O}3.8\ \text{mm}$ Furcation Tubing	\$28.05	Today
FTS3B	Replacement Blade Set (Four Blades) for FTS3	\$8.42	Today
FTS4	Three-Hole Stripping Tool for Fiber Buffers and Jackets	\$29.38	Today
BFG1	Rubber Gripper for Bare Fiber Optics	\$5.00	Today

[Hide Fiber Optic Disposal Unit](#)

Fiber Optic Disposal Unit

[Hide Overview](#)

OVERVIEW

Features

- Safe Solution for Disposing of Cleaved Fiber Ends

- One Circular and One Rectangular Covered Hole in Lid
- Dimensions (L x W x H): 2.75" x 2.75" x 6" (6.99 cm x 6.99 cm x 15.24 cm)
- Disposable: Do Not Reuse

The FTDU offers a safe solution for disposal of cleaved fiber ends, epoxy and syringes, and other sharp objects. This high-density polyethylene (HDPE) bottle features a translucent lid with one circular and one rectangular hole in it (see the photo to the right). Both of these holes have individual covers that snap closed to keep the contents inside the unit. Tabs on the side of the lid hold the covers open for hands-free use. The two hole covers are connected by a plastic tab. Cut, or pull, the covers apart to use the holes separately.



Click to Enlarge

To use, insert sharps into one of the holes in the lid. Do not remove the lid itself from the unit. Use the "Contents" space on the label to clearly indicate what is inside. To protect against spills, keep the holes covered when not inserting sharps. When it is full, securely cover both holes in the lid by inserting the tab on the cover into the slot in the lid and dispose of the canister according to the contents. Do not fill the sharps container past the indicated "Full" line on the label. Do not reuse this unit; it is a one-time-use container only.

Part Number	Description	Price	Availability
FTDU	Fiber Optic Disposal Unit	\$6.95	Today

Visit the *Fiber Optic Termination Tools* page for pricing and availability information:
http://www.thorlabs.com/newgrouppage9_pf.cfm?objectgroup_id=1388