

56 Sparta Avenue • Newton, New Jersey 07860
 (973) 300-3000 Sales • (973) 300-3600 Fax
 www.thorlabs.com



S3FC488 - June 5, 2017

Item # S3FC488 was discontinued on June 5, 2017. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

FIBER-COUPLED LASER SOURCES: VISIBLE, WITH TEC

- ▶ FC/PC Interface for Single Mode Fiber
- ▶ Maximum Output Powers of up to 70 mW
- ▶ Stable Output with Temperature Control



The High-Power Laser Source Front Panel Provides an Enable Button, Display Screen, and Laser Power and Temperature Control



S4FC473
473 nm, 20 mW

[Hide Overview](#)

OVERVIEW

Features

- Two Types of Fabry-Perot Laser Sources
 - Standard Sources for 405 to 520 nm
 - High-Power Sources for 473 to 685 nm with Improved Output Stability and Software Control
- Single Mode, FC/PC Fiber Interface
- Constant Current Operation
- Power Level is Adjustable via Knob and BNC Modulation Input
- Interlock Circuit Provided via 2.5 mm Mono Jack

Thorlabs offers two types of fiber-coupled laser sources with temperature control. The Fabry-Perot Laser Sources with TEC are offered with four different wavelengths, ranging from 405 nm to 520 nm. The High-Power Fabry-Perot Laser Sources with TEC are offered with seven different wavelengths, ranging from 473 nm to 685 nm, and include a software package to control settings and display features. See below for details regarding each type.

The front panel of each FP laser source with TEC includes a display that shows the output power in mW or the temperature in °C, an on/off key, an enable button, a knob to adjust the laser power (drive current), and a remote interlock input (2.5 mm mono jack) for added safety. The standard FP sources have a knob to adjust the temperature (TEC current), while the high-power FP sources have a trimpot to adjust the temperature. The back panel includes a BNC input that allows the laser diode drive current to be controlled via an external voltage source (0 - 5 V) and a remote interlock input. This input enables intensity modulation of the laser source. All of our fiber-pigtailed lasers utilize an angled fiber ferrule at the internal laser/fiber launch point to minimize reflections back

Single Channel Benchtop Laser Sources Selection Guide				
Spectrum	Wavelength	TEC	Laser Type	Fiber Type
Visible	405 - 675 nm	No	Fabry Perot	SM, MM, or PM
	405 - 685 nm	Yes	Fabry Perot	SM
NIR	785 - 1550 nm	No	Fabry Perot	SM or PM
	705 - 1550 nm	Yes	Fabry Perot	SM
	1310 - 1550 nm	Yes	DFB	SM
MIR	2.7 μm	No	Fabry Perot	SM
Other Fiber-Coupled Laser Sources				

into the laser diode, thereby increasing the stability of the laser diode's output. For applications that require several laser sources, consider the temperature-stabilized four-channel fiber coupled laser source.

The high-power FP laser sources permit a maximum modulation frequency up to 100 kHz using a sine wave. A USB connector on the back panel enables computer communication. The included GUI provides control and readout of the laser's output power and setpoints (see *Software* tab).

Please refer to the table to the right for all of our single channel benchtop laser sources.

Key Specifications ^a												
Item #	S3FC405	S3FC473	S3FC488	S3FC520	S4FC473	S4FC488	S4FC520	S4FC637	S4FC642	S4FC660	S4FC685	
Laser Source Type	Standard				High Power							
Wavelength	405 nm	473 nm	488 nm	520 nm	473 nm	488 nm	520 nm	637 nm	642 nm	660 nm	685 nm	
Spectrum												
Full Output Power	Minimum	2.0 mW	5.0 mW	5.0 mW	4.0 mW	18 mW	16 mW	13 mW	65 mW	15 mW	45 mW	13 mW
	Typical	6.0 mW	20.0 mW	20.0 mW	8.0 mW	20 mW	18 mW	15 mW	70 mW	20 mW	50 mW	15 mW
Power Stability	15 Minutes	±0.05 dB	±0.05 dB	±0.05 dB	±0.05 dB	≤0.022 dB (≤0.5%)	≤0.022 dB (≤0.5%)	≤0.009 dB (≤0.2%)	≤0.004 dB (≤0.1%)	≤0.009 dB (≤0.2%)	≤0.004 dB (≤0.1%)	≤0.022 dB (≤0.5%)
	24 Hours	±0.1 dB	±0.1 dB	±0.1 dB	±0.1 dB	≤0.087 dB (≤2.0%)	≤0.087 dB (≤2.0%)	≤0.043 dB (≤1.0%)	≤0.009 dB (≤0.2%)	≤0.009 dB (≤0.2%)	≤0.009 dB (≤0.2%)	≤0.087 dB (≤2.0%)
Fiber Type	Single Mode											

• Complete specifications are available in the *Specs* tab above.

[Hide Specs](#)

S P E C S

FP Laser Source with TEC Specifications					
Item #		S3FC405	S3FC473	S3FC488	S3FC520
Wavelength	Minimum	395 nm	468 nm	483 nm	515 nm
	Typical	405 nm	473 nm	488 nm	520 nm
	Maximum	415 nm	478 nm	488 nm	530 nm
Spectrum ^a					
Full Output Power	Minimum	2.0 mW	5.0 mW	5.0 mW	4.0 mW
	Typical	6.0 mW	20 mW	20 mW	8.0 mW
Laser Class	3B				
Fiber	S405		S460HP		
Output Fiber Connector	FC/PC, Wide 2.1 mm Key Compatible				

• Spectral plots are typical, and actual spectra vary from lot to lot. For further information, please contact Tech Support.

Common Specifications	
Driver	
Power Stability	15 min: ±0.05 dB, 24 hr: ±0.1 dB (After 1 hr Warm-up at 25 ± 10 °C Ambient)
Display Accuracy (mW)	±10% of Actual
Setpoint Resolution	0.01 mW
Adjustment Range	~0 mW to Full Power
TEC	

Stability	0.005 °C / 1 °C
Setpoint Accuracy	±0.25 °C
Setpoint Resolution	±0.1 °C
Adjustment Range	20 ± 1 °C to 30 ± 1 °C
Environmental	
Operating Temperature	15 to 35 °C
Storage Temperature	0 to 50 °C
Input Power	115/230 VAC (Switch Selectable) 50 - 60 Hz
Modulation Input	0 - 5 V = 0 - Full Power, DC or Sinewave Input Only
Modulation Input Impedance	50 Ω
Modulation Bandwidth	5 kHz Full Depth of Modulation 30 kHz Small Signal Modulation

High-Power FP Laser Source with TEC Specifications								
Item #		S4FC473	S4FC488	S4FC520	S4FC637	S4FC642	S4FC660	S4FC685
Wavelength	Minimum	463 nm	478 nm	510 nm	627 nm	632 nm	648 nm	675 nm
	Typical	473 nm	488 nm	520 nm	637 nm	642 nm	658 nm	685 nm
	Maximum	483 nm	498 nm	530 nm	647 nm	652 nm	668 nm	695 nm
Spectrum^a								
Full Output Power^b	Minimum	18 mW	16 mW	13 mW	65 mW	15 mW	45 mW	13 mW
	Typical	20 mW	18 mW	15 mW	70 mW	20 mW	50 mW	15 mW
Laser Class		3B						
Laser Power Stability^c	15 min	≤0.022 dB (≤0.5%)	≤0.022 dB (≤0.5%)	≤0.009 dB (≤0.2%)	≤0.004 dB (≤0.1%)	≤0.009 dB (≤0.2%)	≤0.004 dB (≤0.1%)	≤0.022 dB (≤0.5%)
	24 hr	≤0.087 dB (≤2.0%)	≤0.087 dB (≤2.0%)	≤0.043 dB (≤1.0%)	≤0.009 dB (≤0.2%)	≤0.009 dB (≤0.2%)	≤0.009 dB (≤0.2%)	≤0.087 dB (≤2.0%)
Internal Fiber		460HP			SM600			
Output Fiber Connector		FC/PC, 2.2 mm Wide-Key Slot						

- Spectral plots are typical and actual spectra may vary from lot to lot. Measurements were obtained with an optical spectrum analyzer. For further information, please contact Tech Support.
- Specified at Temperature Set Point of 25 °C
- At 75% of maximum output power. The stability spec was obtained after a 1 hour warm-up period.

Common Specifications		
Driver		
Display Accuracy (mW)	±5% of Actual	
Laser Adjustment Range	Default Setting	Laser Threshold to Max
	Optional Setting	0 mA Current to Max
Laser Adjustment Resolution	0.01 mW (While Laser is Enabled)	
Temperature Adjustment Range	20 °C to 30 °C	
Temperature Setpoint Resolution	0.01 °C	
Modulation Input^a	0 to 5 V (0 to Full Power)	
Modulation Input Impedance	1 kΩ	
Modulation Input Connector^b	BNC	


Modulation Bandwidth^c	≤100 kHz, Sine Wave
General	
Operating Temperature	15 to 35 °C
Storage Temperature	0 to 50 °C
Input Power	85 - 264 VAC 50 - 60 Hz
USB Connector^b	Type B

- Modulation Input voltage directly corresponds to output power, where 5 V = Max Power and 0 V = 0 mW when the front panel knob is set to 0 mW. The maximum voltage will be less than 5 V when using the default laser adjustment range (Laser Threshold to Max).
- The maximum USB and BNC cable length is 3 meters in order to avoid a susceptibility to RF interference according to IEC61000-4-3.
- Waveforms other than sine waves contain components with higher frequency than the overall frequency of the waveform, which may not be followed.


[Hide Temp. Stability](#)

TEMP. STABILITY

The stability data below pertains to our high-power FP laser sources with TEC, and was obtained using an S4FC660 Laser Source. The performance is representative of our other S4FC sources.



A demonstration of the long-term stability of the S4FC660 Laser source, shown over the course of 25 hours. The device is capable of maintaining a constant laser diode temperature, and thus constant output power, over long periods of operation. To obtain the detector signal, the S4FC660 output power was set to 75% of maximum and passed through ND filters reducing the power to 1 mW. The resultant beam was directed into a DET100A detector which was then connected to an oscilloscope via a T-adaptor and high-impedance terminator.




The graph above demonstrates the ability of the S4FC660 to maintain steady temperature regulation even with large changes in the environmental temperature. Even when the environment changes by 20 °C, the S4FC660 is able to maintain a steady and regulated temperature of the laser diode, yielding a consistent power output. To obtain the detector signal, the S4FC660 output power was set to 75% of maximum and passed through ND filters reducing the power to 1 mW. The resultant beam was directed into a DET100A detector which was then connected to an oscilloscope via a T-adaptor and high-impedance terminator.


[Hide Pin Diagrams](#)

PIN DIAGRAMS

Modulation In
BNC Female



Remote Interlock Input
2.5 mm Mono Phono Jack



Terminals must be shorted either by included plug or user device, i.e. external switch, for laser mode "ON" to be

0 to 5 V Max,
50 Ω (S3FC Sources)
1 kΩ (S4FC Sources)

enabled.

USB
USB Type B



Computer Interface
(S4FC Sources)

[Hide Software](#)

SOFTWARE

S4FC Software Package

This software is available only for our S4FC high-power FP laser sources. Our other FP sources cannot be operated by a GUI.

Software

Version 1.1.0


Includes a GUI for control of Thorlabs S4FC Benchtop Laser Sources. To download, click the button below.



GUI Interface

The S4FC software package allows the user to control settings and display features of the benchtop laser source. The settings tab enables several helpful features for controlling the device including setting an output limit, starting the device from the last setting (instead of at 0 mW or threshold power), adjusting the device by current (instead of mW), dimming the display intensity, or enabling the temperature LED to blink when the device is not at its temperature setpoint (making it visually easier to observe when the laser is not at thermal equilibrium).

The controls tab can be used to select the temperature setpoint and power output of the device, as well as enabling the laser's output. Additionally, it will also display the interlock status. When the laser is enabled, this tab will display the current power output and actual temperature of the device.



Click to Enlarge
The Settings tab of the software when the S4FC Benchtop Laser Source is connected.



Click to Enlarge
The Controls tab of the software when the S4FC Benchtop Laser Source is connected and the laser enabled.

[Hide Laser Safety](#)

LASER SAFETY

Laser Safety and Classification

Safe practices and proper usage of safety equipment should be taken into consideration when operating lasers. The eye is susceptible to injury, even from very low levels of laser light. Thorlabs offers a range of laser safety accessories that can be used to reduce the risk of accidents or injuries. Laser emission in the visible and near infrared spectral ranges has the greatest potential for retinal injury, as the cornea and lens are transparent to those wavelengths, and the lens can focus the laser energy onto the retina.

Safe Practices and Light Safety Accessories






- Thorlabs recommends the use of safety eyewear whenever working with laser beams with non-negligible powers (i.e., > Class 1) since metallic tools such as screwdrivers can accidentally redirect a beam.
- Laser goggles designed for specific wavelengths should be clearly available near laser setups to protect the wearer from unintentional laser reflections.
- Goggles are marked with the wavelength range over which protection is afforded and the minimum optical density within that range.
- Laser Safety Curtains, Laser Barriers and Blackout Materials can prevent direct or reflected light from leaving the experimental setup area.
- Thorlabs' Enclosure Systems can be used to contain optical setups to isolate or minimize laser hazards.
- A fiber-pigtailed laser should always be turned off before connecting it to or disconnecting it from another fiber, especially when the laser is at power levels above 10 mW.
- All beams should be terminated at the edge of the table, and laboratory doors should be closed whenever a laser is in use.
- Do not place laser beams at eye level.
- Carry out experiments on an optical table such that all laser beams travel horizontally.
- Remove unnecessary reflective items such as reflective jewelry (e.g., rings, watches, etc.) while working near the beam path.
- Be aware that lenses and other optical devices may reflect a portion of the incident beam from the front or rear surface.
- Operate a laser at the minimum power necessary for any operation.
- If possible, reduce the output power of a laser during alignment procedures.
- Use beam shutters and filters to reduce the beam power.
- Post appropriate warning signs or labels near laser setups or rooms.
- Use laser sign lightboxes if operating Class 3R or 4 lasers (i.e., lasers requiring the use of a safety interlock).
- Do not use Laser Viewing Cards in place of a proper Laser Barrier or Beam Trap.



Laser Classification

Lasers are categorized into different classes according to their ability to cause eye and other damage. The International Electrotechnical Commission (IEC) is a global organization that prepares and publishes international standards for all electrical, electronic, and related technologies. The IEC document 60825-1 outlines the safety of laser products. A description of each class of laser is given below:

Class	Description	Warning Label
1	This class of laser is safe under all conditions of normal use, including use with optical instruments for intrabeam viewing. Lasers in this class do not emit radiation at levels that may cause injury during normal operation, and therefore the maximum permissible exposure (MPE) cannot be exceeded. Class 1 lasers can also include enclosed, high-power lasers where exposure to the radiation is not possible without opening or shutting down the laser.	
1M	Class 1M lasers are safe except when used in conjunction with optical components such as telescopes and microscopes. Lasers belonging to this class emit large-diameter or divergent beams, and the MPE cannot normally be exceeded unless focusing or imaging optics are used to narrow the beam. However, if the beam is refocused, the hazard may be increased and the class may be changed accordingly.	
2	Class 2 lasers, which are limited to 1 mW of visible continuous-wave radiation, are safe because the blink reflex will limit the exposure in the eye to 0.25 seconds. This category only applies to visible radiation (400 - 700 nm).	

2M	Because of the blink reflex, this class of laser is classified as safe as long as the beam is not viewed through optical instruments. This laser class also applies to larger-diameter or diverging laser beams.	
3R	Lasers in this class are considered safe as long as they are handled with restricted beam viewing. The MPE can be exceeded with this class of laser, however, this presents a low risk level to injury. Visible, continuous-wave lasers are limited to 5 mW of output power in this class.	
3B	Class 3B lasers are hazardous to the eye if exposed directly. However, diffuse reflections are not harmful. Safe handling of devices in this class includes wearing protective eyewear where direct viewing of the laser beam may occur. In addition, laser safety signs lightboxes should be used with lasers that require a safety interlock so that the laser cannot be used without the safety light turning on. Class-3B lasers must be equipped with a key switch and a safety interlock.	
4	This class of laser may cause damage to the skin, and also to the eye, even from the viewing of diffuse reflections. These hazards may also apply to indirect or non-specular reflections of the beam, even from apparently matte surfaces. Great care must be taken when handling these lasers. They also represent a fire risk, because they may ignite combustible material. Class 4 lasers must be equipped with a key switch and a safety interlock.	
All class 2 lasers (and higher) must display, in addition to the corresponding sign above, this triangular warning sign		

[Hide Fabry-Perot Laser Sources with Temperature Control](#)

Fabry-Perot Laser Sources with Temperature Control

- ▶ Available Wavelengths: 405 nm, 473 nm, 488 nm, and 520 nm
- ▶ Typical Maximum Output Powers up to 20 mW
- ▶ Stabilities of ±0.05 dB Over 15 Minutes
- ▶ Fiber Patch Cables Sold Separately
- ▶ BNC Connector for Modulating Output with Analog Input
- ▶ Thermoelectric Cooler (TEC) for Temperature Stabilization to within 0.005 °C
- ▶ Custom Wavelengths Available; Contact Tech Support

These Fiber-Coupled Laser Sources have an integrated TEC element that is used to stabilize the temperature of the FP laser diode, which in turn stabilizes the output power and wavelength of the laser diode for a given drive current. The FP laser diode is pigtailed to a single mode fiber that is terminated at an FC/PC bulkhead connector (wide and narrow key compatible) on the front panel.

Also found on the front panel is a display that shows the output power in mW or the temperature in °C, an on/off key, an enable button, a knob to adjust the laser power (drive current), and a knob to adjust the temperature (TEC current). The back panel includes an input that allows the laser diode drive current to be controlled via an external voltage source (0 - 5 V) and a remote interlock input. All of our fiber-pigtailed lasers utilize an angled fiber ferrule at the internal laser/fiber launch point to minimize reflections back into the laser diode, thereby increasing the stability of the laser diode's output. For applications that require several laser sources, consider the temperature-stabilized four-channel fiber coupled laser source.


Note: The laser must be off when connecting or disconnecting fibers from the device, particularly for power levels above 10 mW.

Part Number	Description	Price	Availability
S3FC405	Fabry-Perot Benchtop Laser Source, 405 nm, 2.0 mW (Min), FC/PC	\$2,056.00	Today
S3FC473	Fabry-Perot Benchtop Laser Source, 473 nm, 5.0 mW (Min), FC/PC	\$7,890.00	Lead Time
S3FC520	Fabry-Perot Benchtop Laser Source, 520 nm, 4.0 mW (Min), FC/PC	\$2,260.00	Today

[Hide High-Power Fabry-Perot Laser Sources with Temperature Control](#)

High-Power Fabry-Perot Laser Sources with Temperature Control

- ▶ Available Wavelengths: 473 nm, 488 nm, 520 nm, 637 nm, 642 nm, 660 nm, and 685 nm
- ▶ Typical Maximum Output Powers up to 70 mW
- ▶ Stabilities of ≤0.022 dB Over 15 Minutes at 75% Max Output Power
- ▶ Fiber Patch Cables Sold Separately
- ▶ BNC Connector for Modulating Output with Analog Input
- ▶ High-Current TEC for Temperature Regulation and Stability
- ▶



USB Connector for Computer Control via Included GUI or Command Line

Thorlabs' High-Power Fiber-Coupled Laser Sources are high-performance benchtop Fabry-Perot laser diode sources offering excellent stability and computer control. The internal electronics isolate the laser diode driver from noise coupling (such as that produced by nearby switching supplies). This results in laser sources with lower noise intensity compared to our standard FP Sources with TEC above (see graph to the right). The integrated, high-current TEC element can provide excellent temperature regulation and stability to the laser diode, as illustrated in the graphs on the *Temp. Stability* tab. Even when the ambient temperature changes significantly, the temperature control servo is capable of maintaining constant temperature at the laser diode head. This exceptional temperature stability produces a constant and stable power output from the device.

[Click to Enlarge](#)
 An illustration of the difference in the intensity noise floor of the S3FC520 to the S4FC660. The noise signal from the S4FC660 shows a considerable reduction in the peak-to-peak noise amplitude as well as suppression of high frequency components that are seen in the noise signal of the S3FC520. The S3FC520 plot was vertically offset in order to better compare the two signals.



[Click to Enlarge](#)

The Controls tab of the software when a S4FC Series Benchtop Laser Source is connected and the laser is enabled.

The laser diode is pigtailed to a single mode fiber that is terminated at an FC/PC bulkhead connector (wide and narrow key compatible) on the front panel. An angled fiber ferrule at the internal laser/fiber launch point minimizes reflections back into the laser diode, thereby increasing the stability of the laser diode's output. Also found on the front panel is a display that shows the output power in mW or the temperature in °C, an on/off key, an enable button, a knob to adjust the laser power (drive current), and a trimpot to adjust the temperature (TEC current).

The back panel includes a BNC input that allows the laser diode drive current to be controlled via an external voltage source (0 - 5 V) and a remote interlock input. This input enables intensity modulation of the laser source. Using a sine wave, the output can obtain full-depth modulation at frequencies up to 100 kHz. A USB connector on the back panel enables computer communication. The included GUI provides control and readout of the laser's output power and setpoints (see *Software* tab). The back panel also features a remote interlock input (2.5 mm mono jack) for added safety.

Note that the fiber bulkhead and patch cable ferrule must be cleaned prior to connecting a patch cable. For instructions, please refer to the operating manual. The laser must be off when connecting or disconnecting fibers from the device, particularly for power levels above 10 mW. For applications that require several laser sources, consider the temperature-stabilized four-channel fiber coupled laser source.

Part Number	Description	Price	Availability
S4FC473	Fabry-Perot Benchtop Laser Source, 473 nm, 20 mW (Typ.), FC/PC	\$9,450.00	Today
S4FC488	Fabry-Perot Benchtop Laser Source, 488 nm, 18 mW (Typ.), FC/PC	\$8,720.00	Today
S4FC520	Fabry-Perot Benchtop Laser Source, 520 nm, 15 mW (Typ.), FC/PC	\$3,226.00	Today
S4FC637	Fabry-Perot Benchtop Laser Source, 637 nm, 70 mW (Typ.), FC/PC	\$3,170.00	Today
S4FC642	Fabry-Perot Benchtop Laser Source, 642 nm, 20 mW (Typ.), FC/PC	\$3,170.00	Today
S4FC660	Fabry-Perot Benchtop Laser Source, 660 nm, 50 mW (Typ.), FC/PC	\$2,966.00	Today
S4FC685	Fabry-Perot Benchtop Laser Source, 685 nm, 15 mW (Typ.), FC/PC	\$3,119.00	Today