

RC08FC Off Axis Parabolic Collimator

- Minimize Chromatic Aberrations
- Off-Axis Input Angle: 90°
- Beam Size: Ø8.5 mm
- SM05 External Thread For Easy Mounting
- FC/PC Connector
- Each Collimation Package is Factory Aligned
- Simplifies Free-Space Laser to Fiber Coupling

Specifications

- Clear Aperture for Collimated Beam: Ø11 mm
- Beam Size: Ø8.5 mm
- Numerical Aperture: 0.16
- Effective Focal Length (EFL): 33 mm ±1%
- Parent Focal Length (PFL): 16.5 mm
- Scratch-Dig: 40 - 20
- Reflection (Transmission Through System):
>96% from 400 - 700 nm and 97.5% from 700 - 2000 nm


The RC08FC off-axis collimator creates a 90° angle between the incident and reflected light. The mirror surface is coated with protected silver, which has good reflection properties over a broad wavelength range, making it very useful for systems with many different wavelengths. The RC08FC collimator is pre-aligned to collimate a laser beam propagating from the tip of an FC/PC connectorized fiber with diffraction-limited performance over a large wavelength interval from less than 400 nm to more than 2000 nm. The off-axis parabolic mirror is aligned such that the distance between it and the fiber tip is equal to the focal length of the mirror. The off-axis parabolic collimator is suitable for collimating the output from a single mode fiber. In addition, it can be used to collect light from broadband sources and couple it into a larger core multimode fiber.



RC08FC with Cage Components
See Pages 148 - 154 for 16 mm
Cage System



RC08FC Post Mounted
See Page 238 for
LMR05 Mounted



Please refer to our website for complete models and drawings.

ITEM#	\$	£	€	RMB	DESCRIPTION
RC08FC	\$ 525.00	£ 364.00	€ 466,10	¥ 4,433.20	Ø8.5 mm Beam, Reflective Collimator, FC/PC Connector

Laser Scanning Microscopy and Accessories

Flexible Design for Custom Configurations

Fluorescence Imaging Filter Sets

Dual Wavelength Beam Combiner

Five Position Pinhole Selector

Confocal Base System

ThorVCM Software

Confocal Imaging Controller

4-Channel Laser Source

PMT Unit

Microscopy Stages and Controllers

Objectives and Mounts

See Pages 1394-1405