

SM1 Series of Lens Tube Adapters

- Combine 1/2", 30 mm, and 2" Optics Within the SM1 Series
- Solid One-Piece Construction Offers Excellent Stability

The SM1 Series of Lens Tube Adapters are internally and externally threaded to accept the mating threads of our standard lens tubes. Lens tubes of different sizes can be combined to form optical assemblies.

Adapters with Internal SM1 Threads

ITEM #	\$	£	€	RMB	DESCRIPTION
SM1A1	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41	Adapter with Internal SM1 Threads and External SM05 Threads
SM1A15	\$ 21.00	£ 15.12	€ 18,27	¥ 167.37	Adapter with Internal SM1 Threads and External SM30 Threads
SM2A6	\$ 24.00	£ 17.28	€ 20,88	¥ 191.28	Adapter with Internal SM1 Threads and External SM2 Threads

SM1A6TK

SM1A16

SM1A2

Adapters with External SM1 Threads

ITEM #	\$	£	€	RMB	DESCRIPTION
SM1A6	\$ 18.75	£ 13.50	€ 16,31	¥ 149.44	Adapter with Internal SM05 Threads and External SM1 Threads, 0.15" Thick
SM1A16	\$ 21.00	£ 15.12	€ 18,27	¥ 167.37	Adapter with Internal SM30 Threads and External SM1 Threads
SM1A2	\$ 24.00	£ 17.28	€ 20,88	¥ 191.28	Adapter with Internal SM2 Threads and External SM1 Threads
SM3A1	\$ 29.00	£ 20.88	€ 25,23	¥ 231.13	Adapter with Internal SM3 Threads and External SM1 Threads

Adapter with External SM1 Threads and Increased Thickness

ITEM #	\$	£	€	RMB	DESCRIPTION
SM1A6TK*	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41	Adapter with Internal SM05 Threads and External SM1 Threads, 0.40" Thick

*Two SM05RR Retaining Rings Included

RMS to SM1 Thread Adapters



SM1A3



SM1A4

This mount adapts RMS Threaded Standard microscope objectives to our SM1 series lens tube system.

- SM1A3 features SM1 (1.035"-40) External Threads and RMS (0.800"-36) Internal Threads
- SM1A4 Adapter Designed with SM1 (1.035"-40) Internal Threads and RMS (0.800"-36) External Threads

Low Thermal Conductivity Version (See Page 149)

ITEM #	\$	£	€	RMB
SM1A3	\$ 16.75	£ 12.06	€ 14,57	¥ 133.50
SM1A4	\$ 23.00	£ 16.56	€ 20,01	¥ 183.31

M25 x 0.75 to SM1 Thread Adapters

These adapters are used to integrate M25 x 0.75 objectives into lens tubes or to create an objective with a lens tube.



SM1A11



SM1A12

- SM1A11 Adapter Designed with SM1 (1.035"-40) Internal Threads and M25 x 0.75 External Threads
- SM1A12 has M25 x 0.75 Internal Threads and SM1 (1.035"-40) External Threads

ITEM #	\$	£	€	RMB
SM1A11	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41
SM1A12	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41

C-Mount to SM1 Thread Adapters

These adapters can be used for integrating C-mount and CS-Mount threaded components such as camera lenses into SM1-threaded assemblies. CS-Mount and C-Mount standards are identical except for the flange focal distance.

Low Thermal Conductivity Version (See Page 149)




SM1A9



SM1A10



DC-CS1



ITEM #	\$	£	€	RMB	INTERNAL THREADS	EXTERNAL THREADS
SM1A9	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41	SM1 (1.035"-40)	C-Mount (1.00"-32)
SM1A10	\$ 19.75	£ 14.22	€ 17,18	¥ 157.41	C-Mount (1.00"-32)	SM1 (1.035"-40)
DC-CS1	\$ 26.00	£ 18.72	€ 22,62	¥ 207.22	CS-Mount (1.00"-32)	SM1 (1.035"-40)