

▼ CHAPTERS

Tables/
Breadboards

Mechanics

Optomechanical
Devices

Kits

Lab Supplies

▼ SECTIONS

Optics Supplies

Fiber Supplies

Electrical
Components

Lab Tools

Safety and Blackout

Screws, Nuts,
and Bolts

19" Rack Hardware

Storage


For current pricing,
please see our website.

Guide to Connectorization and Polishing Optical Fibers

- Connectorization
- Cleaving
- Polishing

Free Download at
www.thorlabs.com!

Easy-to-Follow
Step-by-Step
Instructions




ITEM #	\$	£	€	RMB	DESCRIPTION
FN96A	\$ 6.90	£ 4.97	€ 6.00	¥ 54.99	Guide to Connectorization, Polishing, and Cleaving of Fibers

Connectorization Kits

Rubber Polishing Pad (FC and SC Kits only)

- All Connectorization Kits Include
 - Glass Polishing Plate
 - 40 Sheets of Polishing Film
 - Polishing Disc
 - 200X Fiber Scope
 - Diamond Scribe
 - 20 Syringes
 - 2 m Furcation Tubing
 - Epoxy
 - Fiber Stripper
 - Kimwipes®
 - Wash Bottle
 - Crimp Tool
- The CK03 and CK05 Kits Include an Additional Rubber Polishing Pad that is Used to Produce PC (Radiused) Polishes.



Items Not Included

- Connectors (See Pages 412 - 413)
- Fiber (See Pages 1018 - 1064)

ITEM #	\$	£	€	RMB	DESCRIPTION
CK01	\$ 715.12	£ 514.89	€ 622.15	¥ 5,699.51	SMA Connectorization Tool Kit
CK03	\$ 735.42	£ 529.50	€ 639.82	¥ 5,861.30	FC/PC and SC Connectorization Tool Kit
CK05	\$ 771.12	£ 555.21	€ 670.87	¥ 6,145.83	FC/APC Connectorization Tool Kit

Lens Tissues

MC-5 and MC-50E

These extremely soft, premium-grade, 4.9" x 2.9" (124 mm x 73 mm) tissues are shipped in packs of 5 or 50 protective booklets, with 25 sheets per booklet.

ITEM #	\$	£	€	RMB	DESCRIPTION
MC-5	\$ 9.80	£ 7.06	€ 8.53	¥ 78.11	Lens Tissues (5 Booklets)
MC-50E	\$ 78.00	£ 56.16	€ 67.86	¥ 621.66	Lens Tissues (50 Booklets)



Webril Lens Tissues

CP-100

These pure cotton, non-woven, lintless, absorbant cleaning pads are 4" x 8" (100 mm x 200 mm). 100 pads are included per package.

ITEM #	\$	£	€	RMB	DESCRIPTION
CP-100	\$ 11.40	£ 8.21	€ 9.92	¥ 90.86	Webril Pads (100/Pkg)

