

Controller Interface Cables

Thorlabs now offers a comprehensive selection of interface cables that will connect our line of laser diode and TEC controllers to our temperature-controlled laser diode mounts. These high-quality cables are designed for lower noise and higher current applications than standard off-the-shelf cables.



CAB400



CAB420-15



CAB430

Cable and Controller Cross Reference Guide

CONTROLLER SERIES	OEM ITC100	BENCHTOP					PLATFORM			
		LDC200	TED200	LDC340	TED350	ITC500	ITC5000	LDC8000 ¹	TED8000	ITC8000
CAB400		✓		✓		✓		✓		✓
CAB420-15			✓ ³		✓	✓			✓	✓
CAB430	✓									✓ ²

1) Does not apply to LDC8080

2) Applies only to the ITC8000DS15

3) CAB420 is for TED200; CAB420-15 is for TED200C

Our Line of Benchtop Controllers Come Complete With Interface Cables. Use the Cross Reference Guide Above to Determine Replacement Cables.

ITEM#	\$	£	€	RMB	DESCRIPTION
CAB400	\$ 66.00	£ 41.60	€ 61.40	¥ 630.30	Laser Diode Controller Cable
CAB420-15	\$ 72.00	£ 45.40	€ 67.00	¥ 687.60	Temperature Controller Cable
CAB430	\$ 120.00	£ 75.60	€ 111.60	¥ 1,146.00	Y-Cable 9 pin to 15 pin male


Power Cords



T17251



T8003



PC3-5



T2682



T8019

ITEM#	\$ (1-24)*	£ (1-24)*	€ (1-24)*	RMB (1-24)*	DESCRIPTION
T17251	\$ 32.00	£ 20.20	€ 29.80	¥ 305.60	7'6" Power Cord (Package of 5)
T2682	\$ 10.40	£ 6.60	€ 9.70	¥ 99.30	Power Cord with Schuco Plug to IEC
T8003	\$ 26.80	£ 16.90	€ 24.90	¥ 255.90	Power Cord, 2 Way Junction with Schuco Plug to 2 IEC
T8019	\$ 33.60	£ 21.20	€ 31.20	¥ 320.90	Power Cord, 3 Way Junction with Schuco Plug to 3 IEC
PC3-5	\$ 30.00	£ 18.90	€ 27.90	¥ 286.50	3' Power Cord (Package of 5)

*For quantities greater than 24, call for volume discount.

Power Strip

Features

- 6-Outlet Surge Protected Power Strip
- Heavy-Duty Metal Housing
- Illuminated Power Switch
- Surge Protected



SP610



ITEM#	\$	£	€	RMB	DESCRIPTION
SP610	\$ 21.00	£ 13.20	€ 19.50	¥ 200.60	6-Outlet Power Strip

*Call for volume discount.

PDA OEM Power Cable

Features

- Great Solution for OEM Applications
- Black PVC Jacket
- 6-Foot Long Cable
- 22 AWG, 3 Conductors



PDA-C-72

ITEM#	\$	£	€	RMB	DESCRIPTION
PDA-C-72	\$ 16.95	£ 10.70	€ 15.80	¥ 161.90	PDA OEM Power Cable