



Black Poster Board



- 1/16" (1.6 mm) Thick Solid Black Board
- Build Light-Tight Boxes
- Score & Fold - Thick Outer Skin Holds Box Together
- Solid Black Throughout



ITEM#	\$	£	€	RMB	SIZE
TB5	\$ 41.20	£ 26.00	€ 38,30	¥ 393.50	5 Sheets, 20" x 30" (50cm x 76cm)

Black Flocked Self-Adhesive Paper




Image contrast can suffer due to stray light reaching the focal plane of many applications where the light path travels through a tube or other enclosed area. Applying a flat-black paint to the inside may help, but a textured matte-black surface is a much more effective light trap.


Flocked paper is a very black, self-adhesive material made of many fine fibers. When looked at from low angles, it absorbs virtually 100% of the light that strikes it. The fibers will not shed dust or lint, unlike some velvet and felt materials.

- Economical Solution for Eliminating Stray Light

ITEM#	\$*	£*	€*	RMB*	DESCRIPTION
BFP1	\$ 27.80	£ 17.50	€ 25,90	¥ 265.50	Black Flock Paper 30" x 30" Sheet

*Price is for one sheet

Black Masking Tape




General purpose black masking tape is often used along with our blackout materials to form temporary light-tight coverings for sensitive equipment.

ITEM#	\$	£	€	RMB	DESCRIPTION
T137-0.5	\$ 5.89	£ 3.70	€ 5,50	¥ 56.20	1/2" x 60 Yd. (12.7 mm x 55 Meter) Roll
T137-1.0	\$ 7.87	£ 5.00	€ 7,30	¥ 75.20	1.0" x 60 Yd. (25 mm x 55 Meter) Roll
T137-2.0	\$ 13.90	£ 8.80	€ 12,90	¥ 132.70	2.0" x 60 Yd. (50 mm x 55 Meter) Roll

Colored Vinyl Tape

Keep track of lab tools and coordinate application setups by using our colored rolls of vinyl tape to color code fibers, cables, samples, and optomechanics. Besides color coding, the user can make notes on the semi-gloss finish of the vinyl tape with a permanent marker. The temperature range of the tape is from 40°F to 170°F, but it is best to apply the tape at room temperature.



- 1/2" Wide Vinyl Tape
- 36-Yard Roll
- Red, Green, and Yellow Colors

ITEM#	\$	£	€	RMB	DESCRIPTION
VTG-050	\$ 3.50	£ 2.20	€ 3,30	¥ 33.40	Green Vinyl Tape
VTR-050	\$ 3.50	£ 2.20	€ 3,30	¥ 33.40	Red Vinyl Tape
VTY-050	\$ 3.50	£ 2.20	€ 3,30	¥ 33.40	Yellow Vinyl Tape