

56 Sparta Avenue • Newton, New Jersey 07860  
 (973) 300-3000 Sales • (973) 300-3600 Fax  
 www.thorlabs.com

# THORLABS



## MSCL2 - AUG 5, 2020

Item # MSCL2 was discontinued on AUG 5, 2020. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

### MINI-SERIES TABLE CLAMPS

- ▶ Ideal for Securing Assemblies to Mini-Series Breadboard
- ▶ Solid, One-Piece Construction

MSCL3 Table Clamp



MSCL3



MSCL2



MSCL2 Table Clamp

#### OVERVIEW


#### Options

- MSCL3: General-Purpose Clamp (Imperial and Metric Compatible)
- MSCL2(M): Height Adjustable with 4-40 (M3) Cap Screw

Our Mini-Series table clamps are fashioned similarly to our CL series clamps but feature #4 (M3) counterbores. These small, space-conscious clamps are useful when constructing complex optomechanical assemblies on our MS12B(M) breadboard, enabling standard optomechanical components to be secured to the breadboard. The clamps require 4-40 (M3) cap screws and washers, which are conveniently provided in our screw kits.

#### Adjustable Height Clamps

Uniquely, the MSCL2(M) clamp features two cap screws, one for height adjustment and the other for lock-down (see diagram above). This allows for easy adaptation to variable heights without compromising clamping force.



Click to Enlarge

Part Number	Description	Price	Availability
MSCL2/M	Mini-Series Variable Height Clamp, M3 Tapped, Metric	\$8.63	5-8 Days
MSCL3	Mini-Series General-Purpose Clamp	\$8.63	Today
MSCL2	Mini-Series Variable Height Clamp, 4-40 Tapped	\$8.63	5-8 Days

Visit the *Mini-Series Table Clamps* page for pricing and availability information:

[https://www.thorlabs.com/newgroupage9.cfm?objectgroup\\_id=6325](https://www.thorlabs.com/newgroupage9.cfm?objectgroup_id=6325)

[https://www.thorlabs.com/newgroupage9\\_pf.cfm?guide=10&category\\_id=14&objectgroup\\_id=6325](https://www.thorlabs.com/newgroupage9_pf.cfm?guide=10&category_id=14&objectgroup_id=6325)